

**Anti-KLRG-1 antibody [2F1] {FITC} (STJA0004064)**

STJA0004064

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Syrian Hamster monoclonal antibody anti-KLRG-1 is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Syrian Hamster
<b>Reactivity</b>	Human/Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	2F1
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	FITC
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">50928</a>
	<a href="#">10219</a>
<b>Gene Symbol</b>	<a href="#">Klrg1</a>
	<a href="#">KLRG1</a>
<b>Uniprot ID</b>	<a href="#">KLRG1_MOUSE</a>
	<a href="#">KLRG1_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	