

Anti-KLRG-1 antibody [2F1] {PE/Cyanine7} (STJA0004067)

STJA0004067

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Syrian Hamster monoclonal antibody anti-KLRG-1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Syrian Hamster
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2F1
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	50928
	10219
Gene Symbol	Klrg1
	KLRG1
Uniprot ID	KLRG1_MOUSE
	KLRG1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	