

Anti-CD90 antibody [M5/49.4.1] {PE/Cyanine7} (STJA0004070)

STJA0004070

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Thy-1 Membrane Glycoprotein is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M5/49.4.1
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	21838
Gene Symbol	Thy1
Uniprot ID	THY1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	