

Anti-CD19 antibody [1D3] {PE/Cyanine5} (STJA0004071)

STJA0004071

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-B-Lymphocyte antigen Cd19 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1D3
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12478
Gene Symbol	Cd19
Uniprot ID	CD19_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	