

Anti-CD48 antibody [HM48-1] {PerCP} (STJA0004075)

STJA0004075

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Cd48 antigen is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HM48-1
Concentration	5 Mu L/Test
Conjugation	PerCP
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12506
Gene Symbol	Cd48
Uniprot ID	CD48_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	