

Anti-CD54 antibody [YN1/1.7.4] {PerCP/Cyanine5.5} (STJA0004078)

STJA0004078

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Intercellular Adhesion Molecule 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	YN1/1.7.4
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	15894
Gene Symbol	Icam1
Uniprot ID	ICAM1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	