

Anti-CD127/IL-7RA antibody [A7R34] {AF488} (STJA0004080)

STJA0004080

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-7 Receptor Subunit Alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A7R34
Concentration	5 Mu L/Test
Conjugation	AF488
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16197
Gene Symbol	IL7r
Uniprot ID	IL7RA_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	