

Anti-CD272 antibody [PK18.6] {PE} (STJA0004081)

STJA0004081

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-B-And T-Lymphocyte Attenuator is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PK18.6
Concentration	5 Mu L/Test
Conjugation	PE
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	208154
Gene Symbol	Btla
Uniprot ID	BTLA_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	