

Anti-CD28 antibody [37.51] {AF488} (STJA0004082)

STJA0004082

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Syrian Hamster monoclonal antibody anti-T-Cell-Specific Surface Glycoprotein Cd28 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Syrian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	37.51
Concentration	5 Mu L/Test
Conjugation	AF488
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12487
Gene Symbol	Cd28
Uniprot ID	CD28_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	