

## Anti-CD18 antibody [TS1/18.1.2.11] {PE/Cyanine7} (STJA0004118)

STJA0004118

### GENERAL INFORMATION

|                          |                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                          |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Integrin Beta-2 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                          |
| <b>Host/Source</b>       | Mouse                                                                                                       |
| <b>Reactivity</b>        | Human                                                                                                       |

### PRODUCT PROPERTIES

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                                                                   |
| <b>Clone ID</b>            | TS1/18.1.2.11                                                                                                |
| <b>Concentration</b>       | 5 Mu L/Test                                                                                                  |
| <b>Conjugation</b>         | PE/Cyanine7                                                                                                  |
| <b>Purification</b>        |                                                                                                              |
| <b>Dilution Range</b>      |                                                                                                              |
| <b>Formulation</b>         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                  |
| <b>Isotype</b>             | IgG1k                                                                                                        |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |                            |
|---------------------------|----------------------------|
| <b>Gene ID</b>            | <a href="#">3689</a>       |
| <b>Gene Symbol</b>        | <a href="#">ITGB2</a>      |
| <b>Uniprot ID</b>         | <a href="#">ITB2_HUMAN</a> |
| <b>Immunogen</b>          |                            |
| <b>Immunogen Region</b>   |                            |
| <b>Specificity</b>        |                            |
| <b>Immunogen Sequence</b> |                            |