

Anti-CD235 antibody [HIR2] {AF488} (STJA0004126)

STJA0004126

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD235 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HIR2
Concentration	5 µg L/Test
Conjugation	AF488
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	2993
	2994
Gene Symbol	GYPA
	GYPB
Uniprot ID	GLPA_HUMAN
	GLPB_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	