

Anti-CD71 antibody [R17 217.1.3/TIB-219] {PerCP/Cyanine5.5} (STJA0004135)

STJA0004135

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Transferrin Receptor Protein 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R17 217.1.3/TIB-219
Concentration	0.2 mg/mL
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	22042
Gene Symbol	Tfrc
Uniprot ID	TFR1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	