

Anti-CD183/CXCR3 antibody [CXCR3-173] {PerCP/Cyanine5.5} (STJA0004140)

STJA0004140

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-C-X-C Chemokine Receptor Type 3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CXCR3-173
Concentration	0.2 mg/mL
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12766
Gene Symbol	Cxcr3
Uniprot ID	CXCR3_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	