

Anti-CD127/IL-7RA antibody [A019D5] {PE/Cyanine5.5} (STJA0004150)

STJA0004150

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Interleukin-7 Receptor Subunit Alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A019D5
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3575
Gene Symbol	IL7R
Uniprot ID	IL7RA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	