

Anti-CD197/CCR7 antibody [G043H7] {PerCP/Cyanine5.5} (STJA0004151)

STJA0004151

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-C-C Chemokine Receptor Type 7 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G043H7
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	1236
Gene Symbol	CCR7
Uniprot ID	CCR7_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	