

Anti-CD90 antibody [M5/49.4.1] (STJA0004165)

STJA0004165

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Thy-1 Membrane Glycoprotein is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M5/49.4.1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Sterile PBS. < 1.0 EU per mg of the antibody as determined by the LAL method.
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	21838
Gene Symbol	Thy1
Uniprot ID	THY1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	