

Anti-CD34 antibody [581] {PerCP/Cyanine5.5} (STJA0004170)

STJA0004170

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Hematopoietic Progenitor Cell antigen Cd34 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	581
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	947
Gene Symbol	CD34
Uniprot ID	CD34_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	