

Anti-IL-5 antibody [TRFK5] {APC} (STJA0004185)

STJA0004185

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IL-5 is suitable for use in IConfocal MicroscopyC research applications.
Applications	ICFC
Host/Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TRFK5
Concentration	0.2 mg/mL
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3567
	16191
Gene Symbol	IL5
	IL5
Uniprot ID	IL5_HUMAN
	IL5_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	