

Anti-CD15 antibody [HI98] {PerCP} (STJA0004191)

STJA0004191

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD15 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HI98
Concentration	5 Mu L/Test
Conjugation	PerCP
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	2526
Gene Symbol	FUT4
Uniprot ID	FUT4_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	