

Anti-CD83 antibody [HB15e] {PE/TR} (STJA0004203)

STJA0004203

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cd83 antigen is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HB15e
Concentration	5 Mu L/Test
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	9308
Gene Symbol	CD83
Uniprot ID	CD83_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	