

Anti-CD3 antibody [17A2] {PE/TR} (STJA0004213)

STJA0004213

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17A2
Concentration	5 Mu L/Test
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12501 12500 12502 Cd3e Cd3d CD3E_MOUSE CD3D_MOUSE CD3G_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	