

Anti-CD11a antibody [FD441.8] {PE/TR} (STJA0004219)
STJA0004219**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin Alpha-L is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FD441.8
Concentration	0.2 mg/mL
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16408
Gene Symbol	Itgal
Uniprot ID	ITAL_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	