

**Anti-CD32 antibody [IV-3] {PE/TR} (STJA0004239)**

STJA0004239

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Low Affinity Immunoglobulin Gamma Fc Region Receptor II-B is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	IV-3
<b>Concentration</b>	5 Mu L/Test
<b>Conjugation</b>	PE/TR
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG2bk
<b>Storage</b>	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	2213
<b>Gene Symbol</b>	FCGR2B
<b>Uniprot ID</b>	FCG2B_HUMAN
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	