

Anti-CD106 antibody [M/K-2.7] {PE/TR} (STJA0004243)
STJA0004243**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Vascular Cell Adhesion Protein 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M/K-2.7
Concentration	0.2 mg/mL
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	
Gene Symbol	Vcam1
Uniprot ID	VCAM1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	