

Anti-CD45R/B220 antibody [RA3.3A 1/6.1] {PE/TR} (STJA0004251)
STJA0004251**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD45R/B220 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RA3.3A 1/6.1
Concentration	5 Mu L/Test
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	19264 5788
Gene Symbol	Ptpnc PTPRC
Uniprot ID	PTPRC_MOUSE PTPRC_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	