

Anti-CD107a antibody [1D4B] {PerCP/Cyanine5.5} (STJA0004269)

STJA0004269

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Lysosome-Associated Membrane Glycoprotein 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1D4B
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16783
Gene Symbol	Lamp1
Uniprot ID	LAMP1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	