

Anti-CD36 antibody [HM36] {AF647} (STJA0004276)

STJA0004276

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Platelet Glycoprotein 4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HM36
Concentration	5 Mu L/Test
Conjugation	AF647
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12491
Gene Symbol	Cd36
Uniprot ID	CD36_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	