

Anti-CD36 antibody [HM36] {ER780} (STJA0004279)

STJA0004279

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Armenian Hamster monoclonal antibody anti-Platelet Glycoprotein 4 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Armenian Hamster |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | HM36 |
| Concentration | 0.5 mg/ml |
| Conjugation | ER780 |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|------------|
| Gene ID | 12491 |
| Gene Symbol | Cd36 |
| Uniprot ID | CD36_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |