

Anti-Perforin antibody [S16009A] {AF647} (STJA0004285)
STJA0004285**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Perforin-1 is suitable for use in IConfocal MicroscopyC research applications.
Applications	ICFC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S16009A
Concentration	0.5 mg/ml
Conjugation	AF647
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	18646
Gene Symbol	Prf1
Uniprot ID	PERF_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	