

Anti-CD69 antibody [H1.2F3] {PE/Cyanine7} (STJA0004292)

STJA0004292

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Early Activation antigen Cd69 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H1.2F3
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12515
Gene Symbol	Cd69
Uniprot ID	CD69_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	