

Anti-GATD3B antibody (STJA0004294)

STJA0004294

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Putative Glutamine Amidotransferase-Like Class 1 Domain-Containing Protein 3b-Mitochondrial is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications. |
| Applications | WB/IHC/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Antigen affinity purification |
| Dilution Range | WB 1:500-1:2000, IHC 1:30-1:150 |
| Formulation | PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--------------------------------|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | Fusion protein of human GATD3B |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |