

Anti-MFAP3L antibody (STJA0004336)

STJA0004336

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Microfibrillar-Associated Protein 3-Like is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:100-1:300
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9848
Gene Symbol	MFAP3L
Uniprot ID	MFA3L_HUMAN
Immunogen	Recombinant protein of human MFAP3L
Immunogen Region	
Specificity	
Immunogen Sequence	