

Anti-PAPSS2 antibody (STJA0004347)

STJA0004347

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Bifunctional 3 Phosphoadenosine 5 Phosphosulfate Synthase 2 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.02 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9060
Gene Symbol	PAPSS2
Uniprot ID	PAPS2_HUMAN
Immunogen	Fusion protein of human PAPSS2
Immunogen Region	
Specificity	
Immunogen Sequence	