

**Anti-CBX7 antibody (STJA0004358)**

STJA0004358

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Chromobox Protein Homolog 7 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
<b>Applications</b>	WB/IHC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.3 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB 1:500-2000, IHC 1:20-100
<b>Formulation</b>	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">23492</a>
<b>Gene Symbol</b>	<a href="#">CBX7</a>
<b>Uniprot ID</b>	<a href="#">CBX7_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide of human CBX7
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	