

Anti-CELA3A antibody (STJA0004362)

STJA0004362

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Chymotrypsin-Like Elastase Family Member 3a is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.78 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:30-1:150
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10136
Gene Symbol	CELA3A
Uniprot ID	CEL3A_HUMAN
Immunogen	Fusion protein of human CELA3A
Immunogen Region	
Specificity	
Immunogen Sequence	