

Anti-G6PC antibody (STJA0004382)

STJA0004382

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Glucose-6-Phosphatase Catalytic Subunit 1 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.44 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2538
Gene Symbol	G6PC1
Uniprot ID	G6PC1_HUMAN
Immunogen	Synthetic peptide of human G6PC
Immunogen Region	
Specificity	
Immunogen Sequence	