

**Anti-MLLT6 antibody (STJA0004391)**

STJA0004391

**GENERAL INFORMATION**

|                          |                                                                                                                            |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                         |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Protein Af-17 is suitable for use in Immunohistochemistry and ELISA research applications. |
| <b>Applications</b>      | IHC/ELISA                                                                                                                  |
| <b>Host/Source</b>       | Rabbit                                                                                                                     |
| <b>Reactivity</b>        | Human                                                                                                                      |

**PRODUCT PROPERTIES**

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                  |
| <b>Clone ID</b>            |                                                             |
| <b>Concentration</b>       | 1.32 mg/mL                                                  |
| <b>Conjugation</b>         | Unconjugated                                                |
| <b>Purification</b>        | Antigen affinity purification                               |
| <b>Dilution Range</b>      | IHC 1:30-1:150                                              |
| <b>Formulation</b>         | PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4     |
| <b>Isotype</b>             | IgG                                                         |
| <b>Storage Instruction</b> | Store at -20°C for long term, and avoid freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |                               |
|---------------------------|-------------------------------|
| <b>Gene ID</b>            | <a href="#">4302</a>          |
| <b>Gene Symbol</b>        | <a href="#">MLLT6</a>         |
| <b>Uniprot ID</b>         | <a href="#">AF17_HUMAN</a>    |
| <b>Immunogen</b>          | Fusion protein of human MLLT6 |
| <b>Immunogen Region</b>   |                               |
| <b>Specificity</b>        |                               |
| <b>Immunogen Sequence</b> |                               |

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081