

Anti-BPIFB3 antibody (STJA0004395)

STJA0004395

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Bpi Fold-Containing Family B Member 3 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:5000, IHC 1:25-1:100
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	359710
Gene Symbol	BPIFB3
Uniprot ID	BPIB3_HUMAN
Immunogen	Synthetic peptide of human BPIFB3
Immunogen Region	
Specificity	
Immunogen Sequence	