

## Anti-HSPG2 antibody (STJA0004402)

STJA0004402

### GENERAL INFORMATION

|                          |                                                                                                                                                                                   |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                                                                |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Basement Membrane-Specific Heparan Sulfate Proteoglycan Core Protein is suitable for use in Immunohistochemistry and ELISA research applications. |
| <b>Applications</b>      | IHC/ELISA                                                                                                                                                                         |
| <b>Host/Source</b>       | Rabbit                                                                                                                                                                            |
| <b>Reactivity</b>        | Human/Mouse                                                                                                                                                                       |

### PRODUCT PROPERTIES

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                  |
| <b>Clone ID</b>            |                                                             |
| <b>Concentration</b>       | 0.7 mg/mL                                                   |
| <b>Conjugation</b>         | Unconjugated                                                |
| <b>Purification</b>        | Affinity purification                                       |
| <b>Dilution Range</b>      | IHC 1:50-1:200                                              |
| <b>Formulation</b>         | PBS with 0.05% sodium azide and 50% glycerol, PH7.4         |
| <b>Isotype</b>             | IgG                                                         |
| <b>Storage Instruction</b> | Store at -20°C for long term, and avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |                                  |
|---------------------------|----------------------------------|
| <b>Gene ID</b>            | <a href="#">3339</a>             |
| <b>Gene Symbol</b>        | <a href="#">HSPG2</a>            |
| <b>Uniprot ID</b>         | <a href="#">PGBM_HUMAN</a>       |
| <b>Immunogen</b>          | Synthetic peptide of human HSPG2 |
| <b>Immunogen Region</b>   |                                  |
| <b>Specificity</b>        |                                  |
| <b>Immunogen Sequence</b> |                                  |