

Anti-TSPY1 antibody (STJA0004403)

STJA0004403

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Testis-Specific Y-Encoded Protein 1 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.9 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:30-150
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100289087/7258/728137
Gene Symbol	TSPY1
Uniprot ID	TSPY1_HUMAN
Immunogen	Synthetic peptide of human TSPY1
Immunogen Region	
Specificity	
Immunogen Sequence	