

Anti-PLGLB2 antibody (STJA0004404)

STJA0004404

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Plasminogen-Like Protein B is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.14 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5342/5343
Gene Symbol	PLGLB1.PLGLB2
Uniprot ID	PLGB_HUMAN
Immunogen	Fusion protein of human PLGLB2
Immunogen Region	
Specificity	
Immunogen Sequence	