

Anti-AIMP2 antibody (STJA0004416)

STJA0004416

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Aminoacyl Trna Synthase Complex-Interacting Multifunctional Protein 2 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:1000
Formulation	PBS with 0.05% Proclin300 and 50% glycerol, pH7.4.
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7965
Gene Symbol	AIMP2
Uniprot ID	AIMP2_HUMAN
Immunogen	Recombinant Human Aminoacyl tRNA synthase complex-interacting multifunctional protein 2 protein
Immunogen Region	
Specificity	
Immunogen Sequence	