

Anti-PDAP1 antibody (STJA0004418)

STJA0004418

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-28 Kda Heat-And Acid-Stable Phosphoprotein is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.78 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN3 and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11333
Gene Symbol	PDAP1
Uniprot ID	HAP28_HUMAN
Immunogen	Fusion protein of human PDAP1
Immunogen Region	
Specificity	
Immunogen Sequence	