

Anti-ITPR2 antibody (STJA0004427)

STJA0004427

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Inositol 1-4-5-Trisphosphate Receptor Type 2 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.74 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3709
Gene Symbol	ITPR2
Uniprot ID	ITPR2_HUMAN
Immunogen	Synthetic peptide of human ITPR2
Immunogen Region	
Specificity	
Immunogen Sequence	