

## Anti-TNIP1 antibody (STJA0004430)

STJA0004430

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Tnfaip3-Interacting Protein 1 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
<b>Applications</b>	WB/IHC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.8 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB 1:500-2000, IHC 1:25-100
<b>Formulation</b>	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	10318
<b>Gene Symbol</b>	TNIP1
<b>Uniprot ID</b>	TNIP1_HUMAN
<b>Immunogen</b>	Synthetic peptide of human TNIP1
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	