

Anti-PUM3 antibody (STJA0004432)

STJA0004432

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pumilio Homolog 3 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.08 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:1000-1:5000, IHC 1:50-1:200
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9933
Gene Symbol	PUM3
Uniprot ID	PUM3_HUMAN
Immunogen	Synthetic peptide of human PUM3
Immunogen Region	
Specificity	
Immunogen Sequence	