

Anti-ALG11 antibody (STJA0004441)

STJA0004441

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Gdp-Man-Man (3glnac (2-Pp-Dol Alpha-1-2-Mannosyltransferase is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:100-1:300
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	440138
Gene Symbol	ALG11
Uniprot ID	ALG11_HUMAN
Immunogen	Synthetic peptide of human ALG11
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	