

Anti-PLPPR5 antibody (STJA0004445)

STJA0004445

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospholipid Phosphatase-Related Protein Type 5 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.8 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:1000-1:5000, IHC 1:100-1:200
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	163404
Gene Symbol	PLPPR5
Uniprot ID	PLPR5_HUMAN
Immunogen	Synthetic peptide of human PLPPR5
Immunogen Region	
Specificity	
Immunogen Sequence	