

Anti-PLEKHG6 antibody (STJA0004447)

STJA0004447

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pleckstrin Homology Domain-Containing Family G Member 6 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.4 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:100-1:300
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55200
Gene Symbol	PLEKHG6
Uniprot ID	PKHG6_HUMAN
Immunogen	Recombinant protein of human PLEKHG6
Immunogen Region	
Specificity	
Immunogen Sequence	