

Anti-PTPN20 antibody (STJA0004455)

STJA0004455

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tyrosine-Protein Phosphatase Non-Receptor Type 20 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:100-1:300
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26095
Gene Symbol	PTPN20
Uniprot ID	PTN20_HUMAN
Immunogen	Synthetic peptide of human PTPN20A/PTPN20B
Immunogen Region	
Specificity	
Immunogen Sequence	